

Cell viability

How it works

In flow cytometry your data is collected on a per cell basis. So you see information for each cell. This means that a viability dye is a useful addition to exclude dead cells from your data, and when setting up cell sorts.

Viability dyes generally bind to DNA, therefore measure membrane permeability, which is indicative of cell death. However, annexin- V binds to a membrane protein whose abundance increases during apoptosis. Thus, cells which stain positive are dead/undergoing apoptosis.

Reagent list

Annexin V		Protein	Various suppliers
PI	Propidium iodide	DNA	Various
7-AAD	7-amino-actinomycin D	DNA	Various
Hoechst 33342		DNA	Various
FluoroGold		DNA	Molecular probes
LIVE/DEAD viability dyes		Protein	ThermoFisher

Quick info

	Laser	Fixable	Slap it in	Starting point
Annexin V-conjugate	Various	-	No*	Determine
PI	488, 561 nm	-	Yes	2 µg / 10 ⁶ cells
7-AAD	488, 561 nm	-	Yes	10µg / mL
Hoechst 33342	355 nm	Yes	Yes	10µg / mL
FluoroGold	355, 405 nm	-	Yes	10µg / 200 µL
LIVE/DEAD viability dyes	Various	Yes	No*	See kit manual

*kit available

Protocols

PI	7-AAD	Hoechst	FluoroGold
----	-------	---------	------------

- Harvest cells in single suspension (10⁶ cells)
- Wash cells
- Resuspend in stain buffer
- Add stain
- Acquire after 5-10min incubation

Annexin-V	LIVE/DEAD
-----------	-----------

See manufacturers' instructions

Annexin-V: [Thermo](#) | [Merck](#) | [abcam](#)
[LIVE/DEAD](#)